



**Daniela B Husarik**

**Contact**

Daniela B Husarik

## Publications (4)

Winklehner A, Gordic S, Lauk E, Frauenfelder T, Leschka S, Alkadhi H, Husarik D. Automated attenuation-based tube voltage selection for body CTA: Performance evaluation of 192-slice dual-source CT. *Eur Radiol* 2015; 25:2346-53.

Gordic S, Husarik D, Desbiolles L, Leschka S, Frauenfelder T, Alkadhi H. High-pitch coronary CT angiography with third generation dual-source CT: limits of heart rate. *Int J Cardiovasc Imaging* 2014; 30:1173-9.

Husarik D, Saupe N, Pfirrmann C, Jost B, Hodler J, Zanetti M. Ligaments and plicae of the elbow: normal MR imaging variability in 60 asymptomatic subjects. *Radiology* 2010; 257:185-94.

Husarik D, Saupe N, Pfirrmann C, Jost B, Hodler J, Zanetti M. Elbow nerves: MR findings in 60 asymptomatic subjects--normal anatomy, variants, and pitfalls. *Radiology* 2009; 252:148-56.

## Projects (0)

No results found.

---

Kantonsspital St.Gallen

Rorschacher Strasse 95

CH-9007 St.Gallen

T: +41 71 494 11 11

[support.forschung@kssg.ch](mailto:support.forschung@kssg.ch)