



Roger Pfiffner

Contact

Roger Pfiffner

Publications (3)

Abderhalden S, Biro P, Hechelhammer L, Pfiffner R, Pfammatter T. CT-guided navigation of percutaneous hepatic and renal radiofrequency ablation under high-frequency jet ventilation: feasibility study. *J Vasc Interv Radiol* 2011; 22:1275-8.

Hechelhammer L, Rancic Z, Pfiffner R, Mayer D, Meier T, Lachat M, Pfammatter T. Midterm outcome of endovascular repair of ruptured isolated iliac artery aneurysms. *J Vasc Surg* 2010; 52:1159-63.

Pfammatter T, Hechelhammer L, Pfiffner R. A "Buddy Wire" technique for successful OptEase filter retrieval. *J Vasc Interv Radiol* 2009; 20:656-9.

Projects (0)

No results found.

Kantonsspital St.Gallen

Rorschacher Strasse 95

CH-9007 St.Gallen

T: +41 71 494 11 11

support.forschung@kssg.ch